


产品说明书

CD38 Antibody

货号: ABB5745



产品名称	CD38 Antibody
来源宿主	Rabbit
反应种属	Human
克隆类型	Monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB IHC-P/mlHC
使用方法	WB 1:500~1:2000 IHC-P/mlHC 1:100~1:500
有效期	one year
别名	CD38 antigen; ADP ribosyl cyclase 1; Cd38-rs1; ADP ribosyl cyclase; ADP ribosyl cyclase/cyclic ADP-ribose hydrolase; ADP-ribosyl cyclase 1; cADPr hydrolase 1; CD38 antigen (p45); CD38 molecule;
免疫原	A synthesized peptide derived from human CD38
SwissProt	P28907
分子量	45kDa
产品介绍	CD38 is a transmembrane protein involved in several important biological processes, including immune response, insulin secretion, and social behavior. Synthesizes cyclic ADP-ribose, a second messenger for glucose-induced insulin secretion. Also has cADPr

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com