

# 产品说明书

## MUC1 Antibody

货号: ABB5709

产品名称	MUC1 Antibody
来源宿主	Rabbit
反应种属	Human
克隆类型	Monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB IHC-P/mlHC ICC/IF IP FC
使用方法	WB 1:1000~1:2000 IHC-P/mlHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:20 FC 1:20
有效期	one year
别名	MUC1; CA15-3; CD227; DF3 antigen; Episialin; Epithelial membrane antigen; H23 antigen; H23AG; Krebs von den Lungen-6; MAM6; MUC1/ZD; MUC-1/SEC; Mucin-1; Pem; KL-6; MUC-1/X; Tumor-associated mucin;
免疫原	A synthesized peptide derived from human MUC1
SwissProt	P15941
分子量	18-25kDa
产品介绍	This gene is a member of the mucin family and encodes a membrane bound, glycosylated phosphoprotein. This protein is anchored to the apical surface of many epithelia by a transmembrane domain, with the degree of glycosylation varying with cell type. The p

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com