

产品说明书

CD168 Antibody

货号: ABB5688

产品名称	CD168 Antibody
来源宿主	Rabbit
反应种属	Human Mouse Rat
克隆类型	Monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB IHC-P/mlHC IP
使用方法	WB 1:500~1:2000 IHC-P/mlHC 1:50~1:200 IP 1:50
有效期	one year
别名	Hyaluronan mediated motility receptor; Intracellular hyaluronic acid-binding protein; Receptor for hyaluronan-mediated motility; CD168; HMMR; IHABP; RHAMM;
免疫原	A synthesized peptide derived from human CD168
SwissProt	O75330
分子量	85kDa
产品介绍	Involved in cell motility. When hyaluronan binds to HMMR, the phosphorylation of a number of proteins, including PTK2/FAK1 occurs. May also be involved in cellular transformation and metastasis formation, and in regulating extracellular-regulated kinase (

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com