

产品说明书

iNOS Antibody

货号：ABB5993

产品名称	iNOS Antibody
来源宿主	Rabbit
反应种属	Human Mouse
克隆类型	Monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB ICC/IF IP
使用方法	WB 1:500~1:1000 ICC/IF 1:50~1:200 IP 1:50
有效期	one year
别名	HEP-NOS; Hepatocyte NOS; Inducible NO synthase; Inducible NOS; iNOS; NOS2;
免疫原	A synthesized peptide derived from human iNOS
SwissProt	P35228
分子量	131kDa
产品介绍	iNOS Produces nitric oxide (NO) which is a messenger molecule with diverse functions throughout the body. In macrophages, NO mediates tumoricidal and bactericidal actions.

For research use only !

慧蓝生物 竭诚为您服务

地址：上海市浦东新区周浦镇天雄路588弄，电话：19101712317，邮箱：469340997@qq.com