

## 产品说明书

## Cytokeratin 8 Antibody

货号：ABB5383

产品名称	Cytokeratin 8 Antibody
来源宿主	Rabbit
反应种属	Human Mouse
克隆类型	Monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB IHC-P/mlIHC ICC/IF IP FC
使用方法	WB 1:10000~1:20000 IHC-P/mlIHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50
有效期	one year
别名	CARD2; CK-8; CK8; CYK8; cytokeratin 8; Cytokeratin-8; K2C8; K8; keratin 8; Keratin, type II cytoskeletal 8; Keratin-8; KO; KRT8; Type-II keratin Kb8;
免疫原	A synthesized peptide derived from human Cytokeratin 8
SwissProt	P05787
分子量	53kDa
产品介绍	Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle. K8 a type II cytoskeletal keratin. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cell

For research use only!

慧蓝生物 竭诚为您服务

地址：上海市浦东新区周浦镇天雄路588弄，电话：19101712317，邮箱：469340997@qq.com