

产品说明书

N Cadherin Antibody

货号：ABB5015

产品名称	N Cadherin Antibody
来源宿主	Rabbit
反应种属	Human Mouse Rat
克隆类型	Monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB IHC-P/mlHC
使用方法	WB 1:500~1:2000 IHC-P/mlHC 1:50~1:200
有效期	one year
别名	CADH2; CDHN; Cadherin-2; N-cad; N-cadherin; NCAD; Neural-cadherin precursor; cadherin; neural;
免疫原	A synthesized peptide derived from human N Cadherin
SwissProt	P19022
分子量	140kDa
产品介绍	Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. CDH2 may be involved in neuronal rec

For research use only !

慧蓝生物 竭诚为您服务

地址：上海市浦东新区周浦镇天雄路588弄，电话：19101712317，邮箱：469340997@qq.com