

产品说明书

iNOS

货号: ABB00018

产品名称	iNOS
来源宿主	Rabbit
反应种属	Human, Mouse
克隆类型	Rabbit monoclonal
克隆号	03-69
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide. Stable for 12 months.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles. Stable for 12 months.
应用	WB, IHC-P, mlHC, IF-Cell
使用方法	WB: 1:500-1:1000; IHC-P: 1:100-1:200; mlHC: 1:100-1:200; IF-Cell: 1:50-1:100
有效期	one year
别名	NOS2, NOS type II, Inducible NOS, Inducible NO synthase, HEP-NOS
免疫原	Synthetic peptide.
SwissProt	P35228
细胞定位	Cytoskeleton, Nucleus, Cytosol.
分子量	131kDa
产品介绍	NOS2, also named as iNOS and NOS2A, produces nitric oxide (NO) which is a messenger molecule with diverse functions throughout the body. NO is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission, antimicrobial and antitumoral activities. NOS2 is a nitric oxide synthase which is expressed in liver and is inducible by a combination of lipopolysaccharide and certain cytokines.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com