


# 产品说明书

CD163

货号: ABB00024



产品名称	CD163
来源宿主	Rabbit
反应种属	Human Mouse Rat
克隆类型	Rabbit monoclonal
克隆号	R1-64
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	IHC-P, mIHC, WB, IF
使用方法	WB 1:100-1: 1000IHC-P 1:250-1:1000mIHC 1:250-1:1000IF 1:500IP 1:25
有效期	one year
别名	M130, Hemoglobin scavenger receptor
免疫原	N/A
SwissProt	Q86VB7
细胞定位	Secreted, Cell membrane.
分子量	Predicted band size: 125 kDa
产品介绍	CD163 is a monocyte/macrophage specific marker expressed predominantly on cells which possess strong anti-inflammatory potential. The expression of CD163 is strongly induced by anti-inflammatory mediators such as glucocorticoids and interleukin-10, while being inhibited by pro-inflammatory mediators such as interferon-gamma. CD163-expressing mononuclear phagocytes, as well as soluble CD163, may both take part in downregulating an inflammatory response.

For research use only !

**慧蓝生物 竭诚为您服务**

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com