

# 产品说明书

## $\beta$ Actin (Y0519) Rabbit mAb

货号: HL18343

产品名称	$\beta$ Actin (Y0519) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse, Rat,
克隆类型	Rabbit monoclonal
克隆号	Y0519
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 50% Glycerol, 0.05% BSA, 0.05% Proclin 300
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	WB, IHC, mlHC, IF
使用方法	IHC 1:1000-1:10000; mlHC 1:1000-1:10000; WB 1:1000-1:5000; IF 1:100-1:500
有效期	one year
别名	ACTB; Actin; cytoplasmic 1; Beta-actin; Actin $\beta$ ; $\beta$ -Actin; Beta actin; B actin; actin Beta;
SwissProt	P60709
细胞定位	Cytoplasm
分子量	42kD
产品介绍	This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins. [provided by RefSeq, Jul 2008],

For research use only !

**慧蓝生物 竭诚为您服务**

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com