

# 产品说明书

## CD14 (Y0403) Rabbit mAb

货号: HL18248

产品名称	CD14 (Y0403) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse, Rat,
克隆类型	Rabbit monoclonal
克隆号	Y0403
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 50% Glycerol, 0.05% BSA, 0.05% Proclin 300
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	WB, IHC, mIHC, IF
使用方法	IHC 1:1000-1:2500; mIHC 1:1000-1:2500; WB 1:500-1:2500; IF 1:100-1:500
有效期	one year
别名	Monocyte differentiation antigen CD14; Myeloid cell-specific leucine-rich glycoprotein; CD antigen CD14; [Cleaved into: Monocyte differentiation antigen CD14, urinary form; Monocyte differentiation antigen CD14, membrane-bound form]
SwissProt	P08571
细胞定位	Membrane
分子量	50kD
产品介绍	The protein encoded by this gene is a surface antigen that is preferentially expressed on monocytes/macrophages. It cooperates with other proteins to mediate the innate immune response to bacterial lipopolysaccharide. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Mar 2010],

For research use only !

**慧蓝生物 竭诚为您服务**

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com