

# 产品说明书

GRP78/BiP (Y0344) Rabbit mAb

货号: HL18203

产品名称	GRP78/BiP (Y0344) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse, Rat,
克隆类型	Rabbit monoclonal
克隆号	Y0344
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 50% Glycerol, 0.05% BSA, 0.05% Proclin 300
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	WB, IHC, mIHC, IF
使用方法	IHC 1:100-1:500; mIHC 1:100-1:500; WB 1:500-1:2500; IF 1:100-1:500
有效期	one year
别名	78 kDa glucose-regulated protein; GRP-78; Endoplasmic reticulum luminal Ca <sub>2+</sub> -binding protein grp78; Heat shock 70 kDa protein 5; Immunoglobulin heavy chain-binding protein; BiP;
SwissProt	P11021
细胞定位	Cytoplasm
分子量	78kD
产品介绍	The protein encoded by this gene is a member of the heat shock protein 70 (HSP70) family. It is localized in the lumen of the endoplasmic reticulum (ER), and is involved in the folding and assembly of proteins in the ER. As this protein interacts with many ER proteins, it may play a key role in monitoring protein transport through the cell. [provided by RefSeq, Sep 2010],

For research use only!

**慧蓝生物 竭诚为您服务**

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com