

产品说明书

Transferrin Receptor (Y0360) Rabbit mAb

货号: HL18213

产品名称	Transferrin Receptor (Y0360) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse, Rat,
克隆类型	Rabbit monoclonal
克隆号	Y0360
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 50% Glycerol, 0.05% BSA, 0.05% Proclin 300
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	WB, IHC, mIHC, IF
使用方法	IHC 1:100-1:500; mIHC 1:100-1:500; WB 1:500-1:2500; IF 1:100-1:500
有效期	one year
别名	Transferrin receptor protein 1; TR; TfR; TfR1; TfR; T9; p90; CD antigen CD71; [Cleaved into: Transferrin receptor protein 1, serum form; sTfR;]
SwissProt	P02786
细胞定位	Membrane
分子量	84kD
产品介绍	transferrin receptor(TFRC) Homo sapiens This gene encodes a cell surface receptor necessary for cellular iron uptake by the process of receptor-mediated endocytosis. This receptor is required for erythropoiesis and neurologic development. Multiple alternatively spliced variants have been identified. [provided by RefSeq, Sep 2015],

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com