

产品说明书

beta Actin Antibody

货号: ABB0035

产品名称	beta Actin Antibody
来源宿主	Rabbit
反应种属	Human Mouse Rat
克隆类型	Rabbit monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	WB IHC mIHC ICC/IF
使用方法	WB 1:3000~1:10000 IHC 1:50~1:200 ICC/IF 1:50~1:200
有效期	one year
别名	ACTIN; ACTB; BETA-ACTIN; beta actin; actin beta;
免疫原	A synthesized peptide derived from human beta Actin
SwissProt	P60709
参考结果	42kDa
产品介绍	beta-Actin is one of six different actin isoforms that have been identified. The actin molecules found in cells of various species and tissues tend to be very similar in their immunological and physical properties. Therefore, Antibodies against beta-Actin are useful as loading controls for Western Blotting. However it should be noted that levels of beta-Actin may not be stable in certain cells. For example, expression of beta-Actin in adipose tissue is very low and therefore it should not be use

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com