

产品说明书

GAPDH Antibody

货号: ABB0037

产品名称	GAPDH Antibody
来源宿主	Rabbit
反应种属	Human Mouse Rat Chicken fish Monkey
克隆类型	Rabbit monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	WB IHC mIHC ICC/IF IP FC
使用方法	WB 1:3000~1:10000 IHC 1:100~1:500 ICC/IF 1:100~1:250 IP 1:50 FC 1:50
有效期	one year
别名	aging-associated gene 9 protein; G3P; G3PD; GAPD; GAPDH; glyceraldehyde 3-phosphate dehydrogenase; Glyceraldehyde-3-phosphate dehydrogenase; MGC88685
免疫原	A synthesized peptide derived from human GAPDH
SwissProt	P04406
参考结果	36kDa
产品介绍	Glyceraldehyde 3 phosphate dehydrogenase (GAPDH) is well known as one of the key enzymes involved in glycolysis. GAPDH is constitutively abundant expressed in almost cell types at high levels, therefore antibodies against GAPDH are useful as loading controls for Western Blotting. Some pathology factors, such as hypoxia and diabetes, increased or decreased GAPDH expression in certain cell types.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com