

产品说明书

Actin Antibody

货号: ABB0033

产品名称	Actin Antibody
来源宿主	Rabbit
反应种属	Human Mouse Rat Chicken
克隆类型	Rabbit monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	WB IHC mIHC ICC/IF IP FC
使用方法	WB 1:3000~1:10000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50
有效期	one year
别名	a actin; ACTA; ACTA1; Actin alpha skeletal muscle; Actin; actina; actine; aktin; Alpha Actin 1; Alpha skeletal muscle Actin; ASMA; CFTD; CFTD1; CFTDM; MPFD; NEM1; NEM2; NEM3; nemaline myopathy type 3;
免疫原	A synthesized peptide derived from human Actin
SwissProt	P68133
参考结果	42kDa
产品介绍	Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com