

产品说明书

Dylight 800 Goat Anti Rabbit IgG Secondary Antibody

货号: HL19920

产品名称	Dylight 800 Goat Anti Rabbit IgG Secondary Antibody
来源宿主	Goat
反应种属	Rabbit
同种型	IgG(H+L)
标记	Dylight 800
纯化方式	Affinity purified
形式	Liquid
存储溶液	PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.
存储方式	Stable for one year at -15°C to -25°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing
应用	IF, Flow Cyt, WB
使用方法	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications; while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.
有效期	one year
别名	Dylight 800, Goat Anti-Rabbit IgG
产品介绍	This secondary antibodies is Goat Anti-Rabbit IgG (H+L) Antibody Conjugated to DyLight 800, The fluors Excitation Wavelength:777 nm ,Emission Wavelength:794 nm

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com