


产品说明书

pan-CK/pan-Cytokeratin Mouse mAb

货号：HL06815



产品名称	pan-CK/pan-Cytokeratin Mouse mAb
来源宿主	Mouse
反应种属	Human Mouse Rat
克隆类型	Mouse monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	IHC-P, mlIHC
使用方法	IHC-P 1:100-1:400; mlIHC 1:100-1:400
有效期	one year
别名	Cytokeratin Pan, CK pan; CK-pan
免疫原	Synthesized peptide
SwissProt	P13645;P13646;P02533;P19012;P08779;P08727;
细胞定位	Cytoplasmic, Membranous
产品介绍	The antibody can recognize multiple human Cytokeratins, including CK10, 13, 14, 15, 16, 18, 19, and it can be used for immunohistochemical detection of tumors from monolayer and multilayered epithelium.

For research use only!

慧蓝生物 竭诚为您服务

地址：上海市浦东新区周浦镇天雄路588弄，电话：19101712317，邮箱：469340997@qq.com