

产品说明书

HRP conjugated Anti Digoxin antibody

货号: HL19113

产品名称	HRP conjugated Anti Digoxin antibody
反应种属	species independent
标记	HRP conjugated
形式	Liquid
存储溶液	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
存储方式	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
应用	WB, IHC
使用方法	WB 1:500-1:1000; IHC 1:100
产品介绍	Digoxin is a beta blocker drug originally derived from the foxglove plant, Digitalis lanata. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com