

产品说明书

Glut1 (Y0861) Rabbit mAb

货号: HL18630

产品名称	Glut1 (Y0861) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse, Rat
克隆类型	Rabbit monoclonal
克隆号	Y0861
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles. Stable for 12 months.
应用	WB, IHC, mlHC, IF
使用方法	IHC 1:500-1:2500; mlHC 1:500-1:2500; WB 1:5000-1:25000; IF 1:100-1:500
有效期	one year
别名	SLC2A1 ; GLUT1 ; Solute carrier family 2 ; facilitated glucose transporter member 1 ; Glucose transporter type 1, erythrocyte/brain ; GLUT-1 ; HepG2 glucose transporter
SwissProt	P11166;
细胞定位	Cell membrane ; Multi-pass membrane protein . Melanosome .
分子量	Predicted band size: 54 kDa
产品介绍	This gene encodes a major glucose transporter in the mammalian blood-brain barrier. The encoded protein is found primarily in the cell membrane and on the cell surface, where it can also function as a receptor for human T-cell leukemia virus (HTLV) I and II. Mutations in this gene have been found in a family with paroxysmal exertion-induced dyskinesia.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com