

产品说明书

Peroxidase Streptavidin

货号: HL19201

产品名称	Peroxidase Streptavidin
标记	Horseradish Peroxidase (HRP)
形式	Liquid
浓度	1mg/ml
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles. Stable for 12 months.
应用	WB, IHC, ELISA
使用方法	IHC 1:500-1:1000; WB 1:10000-1:100000; ELISA 1:500-1:1000
有效期	one year
产品介绍	Streptavidin, a bacterial protein isolated from Streptomyces avidinii, is similar to egg-white avidin in its ability to bind biotin, and has been used as a replacement for egg-white avidin because of its more favorable chemical properties. Conjugates of streptavidin are recommended for use with Biotin-SP-conjugated antibodies and Biotin-SP-conjugated ChromPure proteins.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com