

产品说明书

Cytokeratin 20

货号: HL06866

产品名称	Cytokeratin 20
来源宿主	Mouse
反应种属	Human, Mouse, Rat
克隆类型	Mouse monoclonal
克隆号	Y044
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles. Stable for 12 months.
应用	IHC, mIHC
使用方法	IHC 1:100-1:200; mIHC 1:100-1:200
有效期	one year
别名	ck20 ;KRT20 ;Keratin, type I cytoskeletal 20 ;Cytokeratin-20 ;CK-20 ;Keratin-20 ;K20 ;Protein IT
免疫原	Synthesized peptide derived from human CK20 AA range: 300-424
SwissProt	P35900
细胞定位	Cytoplasmic, Membranous
产品介绍	The protein encoded by this gene is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This cytokeratin is a major cellular protein of mature enterocytes and goblet cells and is specifically expressed in the gastric and intestinal mucosa. The type I cytokeratin genes are clustered in a region of chromosome 17q12-q21.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com