

产品说明书

HNF-4-alpha

货号: HL26317

产品名称	HNF-4-alpha
来源宿主	Rabbit
反应种属	Human,Mouse,Rat
克隆类型	Rabbit monoclonal
克隆号	H1006
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),0.05% BSA,40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	IHC-P,mIHC,IF-Tissue
使用方法	IHC-P 1:500-1:3000; mIHC 1:500-1:3000; IF-Tissue 1:100-1:200
有效期	one year
别名	HNF4A, HNF4a9, HNF4a8, HNF4a7, HNF4, HNF 4 alpha
免疫原	Recombinant protein within mouse HNF-4-alpha aa 1-250/474.
SwissProt	P41235
细胞定位	Nucleus
产品介绍	The protein encoded by this gene is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com