

产品说明书

Ly6g Rabbit mAb

货号: HL21848

产品名称	Ly6g Rabbit mAb
来源宿主	Rabbit
反应种属	Mouse
克隆类型	Rabbit monoclonal
克隆号	H4001
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Aliquot store at -20°C. Avoid repeated freeze / thaw cycles. Stable for 12 months.
应用	WB; IHC-P; mIHC
使用方法	WB 1: 1000-1: 3000; IHC-P 1:2000-1:5000; mIHC 1:2000-1:5000
有效期	one year
别名	Lymphocyte antigen 6G, Ly-6G, Ly-6G.1, Ly6g
SwissProt	P35461
细胞定位	Cell membrane
分子量	Predicted band size: 14 kDa
产品介绍	Ly6G is a GPI-anchored protein, that is also known as the myeloid differentiation antigen Gr1. The antigen is transiently expressed on monocytes in the bone marrow. The level of antigen expression in the bone marrow directly correlates with granulocyte differentiation and maturation. Ly6G is expressed predominantly on neutrophils, also in a subset of eosinophils, differentiating pre-monocytes, and plasmacytoid dendritic cells.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com