

产品说明书

488 conjugated Goat Anti-Rabbit Secondary Antibody (H+L)

货号: HL19102

产品名称	488 conjugated Goat Anti-Rabbit Secondary Antibody (H+L)
来源宿主	Goat
反应种属	Rabbit
克隆类型	Multiclonal recombinant
标记	488
纯化方式	The antibody was purified from culture media supernatant by immunoaffinity chromatography using Protein G beads.
形式	Liquid
浓度	0.5 mg/mL
存储溶液	PBS with 50% glycerol, 10 mg/mL BSA, 0.1% Proclin-300, pH 7.4.
存储方式	Store at -20°C. Stable for one year after shipment.
应用	IF; FC
使用方法	IF 1:200-1:1000; FC 1:200-1:1000
有效期	one year
产品介绍	This Secondary Antibodies are mixtures of precisely engineered recombinant monoclonal antibodies that recognize multiple complementary epitopes on the same IgG. Each recombinant clone in the final multiclonal mixture is carefully selected after rigorous characterization and screening to ensure the highest level of performance. The selected clones are mixed in various combinations, each of which is then functionally characterized in multiple applications to select the final multiclonal mixture with optimal performance. Amax=493 nm, Emax=522 nm.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com