

产品说明书

PSD95 (Y0455) Rabbit mAb

货号: HL18292

产品名称	PSD95 (Y0455) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse, Rat
克隆类型	Rabbit monoclonal
克隆号	Y0455
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),0.05% BSA,40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	WB, IHC, mlHC, IF
使用方法	IHC 1:250-1:1000; mlHC 1:250-1:1000; WB 1:1000-1:5000; IF 1:100-1:500
有效期	one year
别名	DLG4 ;PSD95 ;Disks large homolog 4 ;Postsynaptic density protein 95 ;PSD-95 ;Synapse-associated protein 90 ;SAP-90 ;SAP90
SwissProt	P78352
细胞定位	Cytoplasm, Membrane
产品介绍	This gene encodes a member of the membrane-associated guanylate kinase (MAGUK) family. It heteromultimerizes with another MAGUK protein, DLG2, and is recruited into NMDA receptor and potassium channel clusters. These two MAGUK proteins may interact at postsynaptic sites to form a multimeric scaffold for the clustering of receptors, ion channels, and associated signaling proteins. Multiple transcript variants encoding different isoforms have been found for this gene.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com