

产品说明书

TRAP (Y1508) Rabbit mAb

货号: HL19350

产品名称	TRAP (Y1508) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse, Rat
克隆类型	Rabbit monoclonal
克隆号	Y1508
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),0.05% BSA,40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	WB, IHC, mIHC, IF
使用方法	IHC 1:100-1:500; mIHC 1:100-1:500; WB 1:100-1:10000; IF 1:50-1:200
有效期	one year
别名	ACP5;Type 5 acid phosphatase;Acid phosphatase 5 tartrate resistant;Acid phosphatase 5, tartrate resistant;ACP5;EC 3.1.3.2;MGC117378;PPA5_HUMAN;serum band 5 tartrate-resistant acid phosphatase;SPENCDI;T5ap;Tartrate resistant acid ATPase;Tartrate resistant acid phosphatase type 5;Tartrate resistant acid phosphatase type 5 precursor;Tartrate-resistant acid ATPase;Tartrate-resistant acid phosphatase type 5;TR AP;TR-AP;TRACP 5;TRAP;TrATPase;Type 5 acid phosphatase
SwissProt	P13686
细胞定位	Cytoplasmic
分子量	Predicted MW: 37 kDa Observed MW: 43 kDa
产品介绍	This gene encodes an iron containing glycoprotein which catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is the most basic of the acid phosphatases and is the only form not inhibited by L(+)-tartrate.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com