

产品说明书

BCL-6 Rabbit mAb(Y1445)

货号: HL19287

产品名称	BCL-6 Rabbit mAb(Y1445)
来源宿主	Rabbit
反应种属	Human, Rat
克隆类型	Rabbit monoclonal
克隆号	Y1445
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),0.05% BSA,40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	WB, IHC, mlHC, IF
使用方法	IHC 1:100-1:500;mlHC 1:100-1:500;WB 1:1000-1:5000; IF 1:100-1:500
有效期	one year
别名	BCL6 ;BCL5 ;LAZ3
SwissProt	P41182;P41183;
细胞定位	Nucleus .
分子量	Predicted MW: 79 kDa Observed MW: 90 kDa
产品介绍	The protein encoded by this gene is a zinc finger transcription factor and contains an N-terminal POZ domain. This protein acts as a sequence-specific repressor of transcription , and has been shown to modulate the transcription of STAT-dependent IL-4 responses of B cells. This protein can interact with a variety of POZ-containing proteins that function as transcription corepressors. This gene is found to be frequently translocated and hypermutated in diffuse large-cell lymphoma (DLCL) , and may be involved in the pathogenesis of DLCL. Alternatively spliced transcript variants encoding different protein isoforms have been found for this gene.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com