

产品说明书

LAMP3 Rabbit mAb(Y1479)

货号: HL19321

产品名称	LAMP3 Rabbit mAb(Y1479)
来源宿主	Rabbit
反应种属	Human
克隆类型	Rabbit monoclonal
克隆号	Y1479
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),0.05% BSA,40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	IHC, mlHC, IF
使用方法	IHC 1:100-1:500;mlHC 1:100-1:500; IF 1:100-1:500
有效期	one year
别名	LAMP3;DCLAMP;TSC403,Lysosome-associated membrane glycoprotein 3 ;Lysosomal-associated membrane protein 3 ;
SwissProt	Q9UQV4;
细胞定位	Cell surface . Lysosome membrane ; Single-pass type I membrane protein . Cytoplasmic vesicle membrane ; Single-pass type I membrane protein . Early endosome membrane ; Single-pass type I membrane protein .
产品介绍	Dendritic cells (DCs) are the most potent antigen-presenting cells. Immature DCs efficiently capture antigens and differentiate into interdigitating dendritic cells (IDCs) in lymphoid tissues that induce primary T-cell responses

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com