

产品说明书

Dkk-1 Rabbit mAb(Y1237)

货号: HL19079

产品名称	Dkk-1 Rabbit mAb(Y1237)
来源宿主	Rabbit
反应种属	Human, Mouse, Rat
克隆类型	Rabbit monoclonal
克隆号	Y1237
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),0.05% BSA,40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.Stable for 12 months.
应用	WB,IHC-P,mIHC,IF
使用方法	IHC 1:100-1:500; mIHC 1:100-1:500; WB 1:1000-1:5000; IF 1:100-1:500
有效期	one year
别名	DKK1 ;Dickkopf-related protein 1 ;Dickkopf-1 ;
SwissProt	O94907; O54908
细胞定位	Secreted
分子量	Predicted MW: 29 kDa Observed MW: 28-40 kDa
产品介绍	This gene encodes a protein that is a member of the dickkopf family. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through its inhibition of the WNT signaling pathway. Elevated levels of DKK1 in bone marrow plasma and peripheral blood is associated with the presence of osteolytic bone lesions in patients with multiple myeloma.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com