

产品说明书

CHAT Rabbit mAb (Y1411)

货号: HL19253

产品名称	CHAT Rabbit mAb (Y1411)
来源宿主	Rabbit
反应种属	Human, Mouse, Rat
克隆类型	Rabbit monoclonal
克隆号	Y1411
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles. Stable for 12 months.
应用	WB, IHC, mlHC, IF, IP, ELISA
使用方法	WB 1:1000-1:5000; IHC 1:100-1:500; mlHC 1:100-1:500; IF 1:100-1:500; ELISA 1:2500-1:10000; IP 1:50-1:200
有效期	one year
别名	CHAT ;Choline O-acetyltransferase ;CHOACTase
SwissProt	P28329;Q03059;P32738;
细胞定位	nucleus ,cytoplasm ,cytosol ,presynapse
分子量	83kDa
产品介绍	This gene encodes an enzyme which catalyzes the biosynthesis of the neurotransmitter acetylcholine. This gene product is a characteristic feature of cholinergic neurons , and changes in these neurons may explain some of the symptoms of Alzheimer's disease. Polymorphisms in this gene have been associated with Alzheimer's disease and mild cognitive impairment. Mutations in this gene are associated with congenital myasthenic syndrome associated with episodic apnea. Multiple transcript variants encoding different isoforms have been found for this gene , and some of these variants have been shown to encode more than one isoform.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com