

产品说明书

CtBP2 Rabbit mAb (Y2140)

货号: HL16429

产品名称	CtBP2 Rabbit mAb (Y2140)
来源宿主	Rabbit
反应种属	Human,Mouse,Rat,
克隆类型	Rabbit monoclonal
克隆号	Y2140
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),50% Glycerol,0.05% BSA,0.05% Proclin 300
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.
应用	WB,IHC,mIHC,IF
使用方法	IHC 1:100-1:500;mIHC 1:100-1:500;WB 1:1000-1:5000;IF 1:100-1:500
有效期	one year
别名	CTBP2 ;C-terminal-binding protein 2
SwissProt	P56545;P56546;Q9EQH5;
细胞定位	Nucleus . Cell junction , synapse .
分子量	Predicted MW: 49 kDa Observed MW: 49 kDa
产品介绍	This gene produces alternative transcripts encoding two distinct proteins. One protein is a transcriptional repressor , while the other isoform is a major component of specialized synapses known as synaptic ribbons. Both proteins contain a NAD ⁺ binding domain similar to NAD ⁺ -dependent 2-hydroxyacid dehydrogenases. A portion of the 3' untranslated region was used to map this gene to chromosome 21q21.3; however , it was noted that similar loci elsewhere in the genome are likely. Blast analysis shows that this gene is present on chromosome 10. Several transcript variants encoding two different isoforms have been found for this gene.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com