

产品说明书

CD27 Rabbit mAb (Y1331)

货号: HL19173

产品名称	CD27 Rabbit mAb (Y1331)
来源宿主	Rabbit
反应种属	Human,Mouse,Rat,
克隆类型	Rabbit monoclonal
克隆号	Y1331
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),50% Glycerol,0.05% BSA,0.05% Proclin 300
存储方式	Store at +4℃ after thawing. Aliquot store at -20℃. Avoid repeated freeze / thaw cycles.
应用	WB,IHC,mIHC,IF
使用方法	IHC 1:100-1:500;mIHC 1:100-1:500;WB 1:5000-1:25000;IF 1:100-1:500
有效期	one year
别名	CD27 TNFRSF7
SwissProt	P26842;P41272;
细胞定位	Membrane; Single-pass type I membrane protein.
分子量	Predicted MW: 29 kDa Observed MW: 55 kDa
产品介绍	The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is required for generation and long-term maintenance of T cell immunity. It binds to ligand CD70 , and plays a key role in regulating B-cell activation and immunoglobulin synthesis. This receptor transduces signals that lead to the activation of NF-kappaB and MAPK8/JNK. Adaptor proteins TRAF2 and TRAF5 have been shown to mediate the signaling process of this receptor. CD27-binding protein (SIVA) , a proapoptotic protein , can bind to this receptor and is thought to play an important role in the apoptosis induced by this receptor.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com